

Tony: a tool for melody transcription - Feature #903 integrate and test Constant-Q spectrogram

2014-03-12 10:20 AM - Chris Cannam

Status:	New	Start date:	2014-03-12
Priority:	Normal	Due date:	
Assignee:	Chris Cannam	% Done:	0%
Category:		Estimated time:	0.00 hour
Target version:			
Description			
<p>Integrate a Constant-Q spectrogram and compare it with the current spectrogram.</p> <p>This is dependent on sorting out the y-axis scaling in the Colour 3D Plot layer so that it supports bins with non-linear value increments as an auto-align target for other layers.</p>			